



# Biomedical Imaging & Analysis

Lecture 5. Fall 2014

*Image Formation & Visualization (V):*  
***Optical Imaging Systems.***


*Prahlad G Menon, PhD*

*Assistant Professor*

*Sun Yat-sen University – Carnegie Mellon University (SYSU-CMU)*

*Joint Institute of Engineering*

# Basic Microscope Setup




- How a Lens Works
- Diffraction Pattern
- Linear System's view of an optical system

# Performance Metrics of a Microscope

- Resolution: the smallest distance that can be resolved.
- Field of view: the area of a specimen that can be observed and recorded in an image.
- Depth-of-field: the axial distance (depth) in the specimen that appears in focus in an image.
- Light collection power: it determines image brightness.
  - Image bit depth: A 16 bit images have a finer discretization of the range of intensity data stored than an 8 bit image, at each pixel / voxel. Without recording ability the light collection power of a microscope (or telescope for that matter) cannot be taken advantage of...

# Optics: Some Basics (I)

- Lens and Magnification (courtesy, Wikipedia)
- Thin lens approximation (a simple math model)
- Real v/s Imaginary Images




# Lenses of the Compound Microscope

Gives the viewer an erect enlarged virtual image.

Oculars (eyepieces): provide some magnification, focus image at eye (can be 1 or 2 of them)

Objectives: located on rotating nosepiece, provide magnification and resolving power (may be 3 to 6 of them)

Condenser: located under stage, focuses light on the specimen. It may also contain an iris diaphragm that controls the size of the cone of light entering the condenser.



# Magnification

The overall magnification is given as the product of the lenses and the distance over which the image is projected:


$$M = \frac{D \cdot M_1 \cdot M_2}{250mm}$$

where:

D = projection (tube) length (usually = 250 mm);

M<sub>1</sub>, M<sub>2</sub> = magnification of objective and ocular.

250 mm = minimum distance of distinct vision for 20/20 eyes.



# Photon Energy (II)

- Energy of a photon: Planck's law

$$E = h\nu = h \frac{c}{\lambda}$$

$h$ : Planck's constant;  $6.626 \times 10^{-34} \text{ J}\cdot\text{s}$


$\nu$ : frequency of light;

$\lambda$ : wavelength of light

$c$ : speed of light

energy of a photon =  
 $3.973 \times 10^{-19} \text{ J}$  at 500nm

- Shorter waves have higher energy.



# Microscope Modeling – breaking optical imaging down to a science...

What is a Linear System..?  $\mathbf{Ax} = \mathbf{b}$

*Linear transforms you probably know :*

- Rotation, translation, scaling – spatial transforms
- Convolution

*Wolfram – definition of a linear operator:*

A linear transformation between two vector spaces  $V$  and  $W$  is a map  $T : V \rightarrow W$  such that

1.  $T(\mathbf{v}_1 + \mathbf{v}_2) = T(\mathbf{v}_1) + T(\mathbf{v}_2)$  for any vectors  $\mathbf{v}_1$  and  $\mathbf{v}_2$  in  $V$ , and
2.  $T(\alpha \mathbf{v}) = \alpha T(\mathbf{v})$  for any scalar  $\alpha$ .



# A Microscope can be modeled as a Linear System








Figure 4

**Today's  
class**




# Review: Basic Concepts

- An image records spatiotemporal information of a biological process.
- An image can be considered as both a matrix and a surface.



# Recording the Image: what is actually happening at each Pixel?



# Camera Noise Model

- Signal  $S = I \cdot QE \cdot T$

$T$ : exposure/integration time  
 $QE$ : quantum efficiency

- Signal shot noise

$$N_{shot} = \sqrt{S}$$


- Camera noise

$$N_{dark} = \sqrt{D \cdot T} \quad N_{camera} = \sqrt{N_{read}^2 + N_{dark}^2}$$

$D$ : dark current

- Total noise

$$N_{total} = \sqrt{N_{shot}^2 + N_{read}^2 + N_{dark}^2}$$



# Read Noise

- Minimized by careful electronic design
- Under low-light/low-signal conditions where read noise exceeds photon noise, data is read noise limited
- Read noise is not as significant in high-signal applications
- Read noise =  $\text{std} \times \text{system gain} \times 0.707$

(std of subtracted bias images)



Reading all the buckets -  
what's my Error?

# Dark Current

---

- Dark Current is created by heat and cosmic noise and can be reduced by cooling
- Dark Current builds over time unlike read noise
- Dark current reduction is sensor dependent
- For example, some sensors will halve dark current for every 7 degrees of cooling; some require more cooling
- Other technologies can be applied which reduce the cooling required

*Retiga SRV (cooled to -30) Dark Current 0.15 e/p/s*

*Exi Blue (cooled to zero) Dark Current 0.005 e/p/s*

---